


A hook and loop strip is supplied with each Neo-Strip enables attachment of the cannula



The hook and loop is wrapped around the cannula in order to secure to Neo-Strip



The Neo-Strip is made from Hydrocolloid material with an adhesive backing for gently attaching to the baby's delicate skin, without causing damage, but of sufficient retention to allow insertion of cannula or feeding tubes

- Developed to replace adhesives that can be harmful to the babies delicate skin
- The Neo-Strip system offers repositioning through the hook and loop strip
- The hook and loop strip can be wrapped around the tube and attached to the Neo-Strip to secure the cannula
- Hydrocolloid material gentle to skin

## SIZES AVAILABLE

When ordering please quote the following (minimum qty 1 box):

DESCRIPTION	ORDER#	QTY/BOX
Neo-Strip	NS001	50

**Conlett**

U25/17 Cairns St  
Loganholme  
QLD 4129  
Australia

Tel: (+61) 7 3209 9957  
Email: [contact@conlett.com](mailto:contact@conlett.com)